



An ideal supplement for your MOTomed trainer

A good trainer is essential but so is the right training chair. A VELA exercise chair gives you the best conditions to exercise or retrain your arms and legs safely and ergonomically.

- ✓ Train safely and ergonomically
- ✓ Easy to seat, transport and position the user
- ✓ Easy to adjust to the individual user's needs

Train safely and ergonomically

A VELA exercise chair is a safe training chair due to its wheel brake and anti-tip bar. These safety features provide extra stability so the chair cannot tip backwards even with powerful movement. It is an ergonomic chair with adjustable backrest and comfortable padding that provide optimum training conditions for the user.

Easy to seat, transport and position the user

Patients can easily be seated, transported and positioned at your MOTomed trainer due to the free-rolling castors, which are securely locked after positioning.

Easy to adjust to the individual user's needs

The chair is easy to adjust for the individual user to provide the right sitting height and back support. The seat tilt is an efficient way to get the right hip angle for optimum movement of the legs.





VELA 'Medium' Exercise Chair

Item no. 80150

The VELA 'Medium' Exercise Chair has a low backrest that fully supports the lower back while allowing free movement of the back, shoulders and arms during training.

Seat height:	40 - 55 cm
<i>Can be adjusted in the following intervals</i>	42 - 57 cm
	50.5 - 65.5 cm
Seat (W x D):	44 x 44 cm
Backrest (W x H):	37 x 30 cm
Frame size:	55 x 55 cm
Max user weight:	160 kg



VELA 'Large' Exercise Chair

Item no. 80151

The VELA 'Large' Exercise Chair has a higher backrest and larger seat than the 'Medium' model. In addition, the backrest can be locked at a specific angle.

Seat height:	40 - 55 cm
<i>Can be adjusted in the following intervals</i>	42 - 57 cm
	50.5 - 65.5 cm
Seat (W x D):	47 x 49 cm
Backrest (W x H):	43 x 45 cm
Frame size:	55 x 55 cm
Max user weight:	160 kg

FABRICS AND COLOURS

The VELA exercise chairs are upholstered in compound textile by Maharam Design Studio. Compound provides a comfortable surface that does not get slippery or sweaty even when sitting down for a longer time.

We have the following colours on stock but can provide all colours if needed at an additional cost.



CLEANING

Wipe regularly with warm water and neutral detergent. Use a clean white cloth or sponge.

DISINFECTION

Use detergents containing up to 70-85% ethanol or chlorine disinfection containing active chlorine at 1000-1200 ppm. After cleaning, rinse repeatedly with cold, clean water.

Danish Quality

VELA is a Danish family-owned company founded in 1935 by Vermund Larsen. VELA is headquartered in Aalborg, Denmark and has its own development department, and production and padding workshop. In order to be close to our customers, we have a strong, international dealer network.

CHAIRS DESIGNED WITH CARE AND PRECISION

The development department at VELA is constantly focused on optimising safety, functionality and ergonomics in cooperation with health professionals. We take into consideration that chairs in the healthcare sector are exposed to severe wear and tear for many hours each day. At VELA, we use hard-wearing mechanical parts and extra-durable seating materials to keep the chairs clean and functional for many years. All VELA chairs are CE marked and made in Denmark.



Dealer information